

KERN 337-04

Individual weight, OIML Class F2, 10 g, Knob, stainless steelfine turned (OIML)



Services (optional)

| | |
|--|---------|
| Article number for DAkkS calibration | 962-434 |
| Article number for guaranteed plus-deviation | 921-434 |

Pictograms

STANDARD



Category

| | |
|------------------|-------------------|
| Brand | KERN |
| Product category | Individual weight |
| Product group | OIML F2 |
| Product family | 337 |

Measuring System

| | |
|---|--------|
| Conformity assessment possible | no |
| Test weight - nominal value | 10 g |
| Test weight - volume (cm ³) | 1,25 |
| Test weight - volume uncertainty (cm ³) | 0,0156 |
| Permissible tolerance according to OIML (+/-) | 0,6 mg |
| Test weight - Density [Min] (OIML) | 4000 |
| Weight OIML class | F2 |

Construction

| | |
|------------------|-----------------------------------|
| Shape weight | Knob |
| Dimensions (Ø×H) | 10×19 mm |
| Material | stainless steelfine turned (OIML) |

Packing & Shipping

| | |
|------------------------------|----------------|
| Readability force [d] (N) | 1 d |
| Dimensions packaging (W×D×H) | 10×10×19 mm |
| Net weight | 0,01 kg |
| Shipping method | Parcel service |
| Net weight approx. | 0,05 kg |
| Gross weight approx. | 0,05 kg |
| Shipping weight | 0,1 kg |