

KERN 337-05

Individual weight, OIML Class F2, 20 g, Knob, stainless steelfine turned



Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	13×13×23 mm
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Optional cases

Case - Plastic article number	347-050-400
Case - Aluminium article number	317-050-600
Case - Wood article number	337-050-200

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F2
Product family	337

Measuring System

Conformity assessment possible	no
Test weight - volume (cm ³)	2,5
Test weight - volume uncertainty (cm ³)	0,0313
Test weight - Density (kg/m ³)	7950
Test weight - Density uncertainty (kg/m ³)	140
Test weight - Density [Min] (OIML)	4800
Test weight - Density [Max] (OIML)	24000
Test weight - nominal value	20 g
OIML-compliant	✓
Weight OIML class	F2
Permissible tolerance according to OIML (+/-)	0,8 mg

Construction

Shape weight	Knob
Dimensions (Ø×H)	13×23,5 mm
Material	stainless steelfine turned

Services (optional)

Article number for DAkkS calibration	962-435
Article number for guaranteed plus-deviation	921-435