

KERN 337-080-200

KERN

Wooden box for individual weights (Cylindrical/compact, 200 g)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for OIML classes	F2 + M1
Case suitable for nominal weight values	200 g
Suitable weight shape for case	Cylindrical/compact

Construction

Dimension (W×D×H)	60×60×66 mm
Material	beech
Case - Smallest possible weight	200 g
Case - Highest possible weight	200 g

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	60×60×66 mm
Net weight	0,19 kg
Shipping method	Parcel service
Net weight approx.	0,20 kg
Gross weight approx.	0,20 kg
Shipping weight	0,19 kg

Pictograms

STANDARD

