

Product brochure

07.08.2026

KERN 337-11

Individual weight, OIML Class F2, 1 kg, Knob, stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
315-245	Tweezers 315-245 (stainless steel 250 mm)
317-110-600	Aluminium weight case 317-110-600 (for nominal values 1 kg for classes E1 - M3 for design Cylindrical/compact)
317-280	Gloves 317-280 (cotton, 1 pair)
337-110-200	Wood weight case 337-110-200 (for nominal values 1 kg for classes F2 + M1 for design Cylindrical/compact)
347-110-400	Plastic weight case 347-110-400 (for nominal values 1 kg for classes F1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
921-441	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-441
962-441	DAkKS calibration certificate 962-441
962-441R	DAkKS Recalibration 962-441R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days