

Product brochure

07.08.2026

KERN 337-12

Individual weight, OIML Class F2, 2 kg, Knob, stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
315-245	Tweezers 315-245 (stainless steel 250 mm)
317-120-600	Aluminium weight case 317-120-600 (for nominal values 2 kg for classes E1 - M3 for design Cylindrical/compact)
317-280	Gloves 317-280 (cotton, 1 pair)
337-120-200	Wood weight case 337-120-200 (for nominal values 2 kg for classes F2 + M1 for design Cylindrical/compact)
347-120-400	Plastic weight case 347-120-400 (for nominal values 2 kg for classes F1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
921-442	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-442
962-442	DAkkS calibration certificate 962-442
962-442R	DAkkS Recalibration 962-442R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days