

Product brochure

07.08.2026

KERN 337-120-200

Wooden box for individual weights (Cylindrical/compact, 2 kg)

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



FEATURES

- Suitable for individual weights KERN 337, 347, 357, 367

PICTOGRAMS

Standard



SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day