

Product brochure

07.08.2026

KERN 337-13

Individual weight, OIML Class F2, 5 kg, Knob, stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

PICTOGRAMS

Standard



Option



ACCESSORIES

Model	Description
317-130-600	Aluminium weight case 317-130-600 (for nominal values 5 kg for classes E1 - M3 for design Cylindrical/compact)
317-280	Gloves 317-280 (cotton, 1 pair)
317-290	Gloves 317-290 (leather, cotton, 1 pair)
337-130-200	Wood weight case 337-130-200 (for nominal values 5 kg for classes F2 + M1 for design Cylindrical/compact)
347-130-400	Plastic weight case 347-130-400 (for nominal values 5 kg for classes F1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
921-443	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-443
962-443	DAkks calibration certificate 962-443
962-443R	DAkks Recalibration 962-443R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days