

KERN 337-13



Individual weight, OIML Class F2, 5 kg, Knob, stainless steelfine turned



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F2
Product family	337

Measuring System

Conformity assessment possible	no
Test weight - volume (cm ³)	625
Test weight - volume uncertainty (cm ³)	7,8125
Test weight - Density (kg/m ³)	7950
Test weight - Density uncertainty (kg/m ³)	140
Test weight - Density [Min] (OIML)	6400
Test weight - Density [Max] (OIML)	10700
Test weight - nominal value	5 kg
OIML-compliant	✓
Weight OIML class	F2
Permissible tolerance according to OIML (+/-)	80 mg

Construction

Shape weight	Knob
Dimensions (Ø×H)	80×147 mm
Material	stainless steelfine turned

Services (optional)

Article number for DAkkS calibration	962-443
Article number for guaranteed plus-deviation	921-443

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	80×80×147 mm
Shipping method	Parcel service
Net weight approx.	5,0 kg
Gross weight approx.	5,0 kg
Shipping weight	5 kg

Optional cases

Case - Plastic article number	347-130-400
Case - Aluminium article number	317-130-600
Case - Wood article number	337-130-200

Pictograms

STANDARD

