

# KERN 337-140-200

**KERN**

Wooden box for individual weights (Cylindrical/compact, 10 kg)



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

## Measuring System

Case suitable for OIML classes	F2 + M1
Case suitable for nominal weight values	10 kg
Suitable weight shape for case	Cylindrical/compact

## Construction

Dimension (W×D×H)	126×126×200 mm
Material	beech
Case - Smallest possible weight	10 kg
Case - Highest possible weight	10 kg

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	126×126×200 mm
Net weight	0,82 kg
Shipping method	Parcel service
Net weight approx.	0,85 kg
Gross weight approx.	0,85 kg
Shipping weight	0,82 kg

## Pictograms

STANDARD

