

KERN 337-141-200

KERN

Wooden box for individual weights (Check weight, 10 kg)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for OIML classes	F2 + M1
Case suitable for nominal weight values	10 kg
Suitable weight shape for case	Check weight

Construction

Dimension (W×D×H)	205×205×190 mm
Material	beech
Case - Smallest possible weight	10 kg
Case - Highest possible weight	10 kg

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	205×205×190 mm
Net weight	3,1 kg
Shipping method	Parcel service
Net weight approx.	3,2 kg
Gross weight approx.	3,2 kg
Shipping weight	3,1 kg

Pictograms

STANDARD

