

KERN 337-150-200

KERN

Wooden box for individual weights (Cylindrical/compact, 20 kg)



Category

| | |
|------------------|---------------|
| Brand | KERN |
| Product category | Cases & Boxes |
| Product group | Wooden box |

Measuring System

| | |
|---|---------------------|
| Case suitable for OIML classes | F2 + M1 |
| Case suitable for nominal weight values | 20 kg |
| Suitable weight shape for case | Cylindrical/compact |

Construction

| | |
|---------------------------------|----------------|
| Dimension (W×D×H) | 185×185×300 mm |
| Material | beech |
| Case - Smallest possible weight | 20 kg |
| Case - Highest possible weight | 20 kg |

Packing & Shipping

| | |
|------------------------------|----------------|
| Delivery time | 1 d |
| Dimensions packaging (W×D×H) | 185×185×300 mm |
| Net weight | 5,2 kg |
| Shipping method | Parcel service |
| Net weight approx. | 6 kg |
| Gross weight approx. | 6 kg |
| Shipping weight | 5,2 kg |

Pictograms

STANDARD

