

KERN 337-151-100

KERN

Wooden box for individual weights (Check weight, 20 kg)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for OIML classes	F1
Case suitable for nominal weight values	20 kg
Suitable weight shape for case	Check weight

Construction

Dimension (W×D×H)	205×205×290 mm
Material	beech
Case - Smallest possible weight	20 kg
Case - Highest possible weight	20 kg

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	205×205×290 mm
Net weight	4 kg
Shipping method	Parcel service
Net weight approx.	4,0 kg
Gross weight approx.	4,0 kg
Shipping weight	4 kg

Pictograms

STANDARD

