

# KERN 337-151-200

**KERN**

Wooden box for individual weights (Check weight, 20 kg)



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

## Measuring System

Case suitable for OIML classes	F2 + M1
Case suitable for nominal weight values	20 kg
Suitable weight shape for case	Check weight

## Construction

Dimension (W×D×H)	205×205×290 mm
Material	beech
Case - Smallest possible weight	20 kg
Case - Highest possible weight	20 kg

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	205×205×290 mm
Net weight	4 kg
Shipping method	Parcel service
Net weight approx.	4,0 kg
Gross weight approx.	4,0 kg
Shipping weight	4 kg

## Pictograms

STANDARD

