

Product brochure

07.08.2026

# KERN 337-160-200

Wooden box for individual weights (Cylindrical/compact, 50 kg)

ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

+49 7433 9933-0  
info@kern-sohn.com  
www.kern-sohn.com

FEATURES

- Suitable for individual weights KERN 337, 347, 357, 367

PICTOGRAMS

Standard



ACCESSORIES

| Model   | Description                           |
|---------|---------------------------------------|
| 970-102 | Warranty extension (+2 years) 970-102 |

SERVICES

| Model   | Description                                                  |
|---------|--------------------------------------------------------------|
| 920-350 | Marking by laser engraving (individual marking) KERN 920-350 |

ENGLISH

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day