

KERN 337-160-200

KERN

Wooden box for individual weights (Cylindrical/compact, 50 kg)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for OIML classes	F2 + M1
Case suitable for nominal weight values	50 kg
Suitable weight shape for case	Cylindrical/compact

Construction

Material	beech
Case - Smallest possible weight	50 kg
Case - Highest possible weight	50 kg

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	10×10×10 mm
Net weight	7,2 kg
Shipping method	Parcel service
Net weight approx.	8 kg
Gross weight approx.	8 kg
Shipping weight	7,2 kg

Pictograms

STANDARD

