

KERN 337-161-100

KERN

Wooden box for individual weights (Check weight, 50 kg)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for OIML classes	F1
Case suitable for nominal weight values	50 kg
Suitable weight shape for case	Check weight

Construction

Dimension (W×D×H)	290×290×325 mm
Material	beech
Case - Smallest possible weight	50 kg
Case - Highest possible weight	50 kg

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	290×290×325 mm
Net weight	7 kg
Shipping method	Parcel service
Net weight approx.	7 kg
Gross weight approx.	7 kg
Shipping weight	7 kg

Pictograms

STANDARD

