

KERN 338-090-200

KERN

Wooden box for individual weights (Cylindrical/compact, 1 mg - 500 mg (einzel))



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for OIML classes	E1 - M2
Case suitable for nominal weight values	1 mg - 500 mg (einzel)
Suitable weight shape for case	Cylindrical/compact

Construction

Dimension (W×D×H)	60×60×35 mm
Material	beech
Case - Smallest possible weight	1 mg
Case - Highest possible weight	500 mg

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	60×60×35 mm
Net weight	0,13 kg
Shipping method	Parcel service
Net weight approx.	0,15 kg
Gross weight approx.	0,15 kg
Shipping weight	0,13 kg

Pictograms

STANDARD

