

07.08.2026

Weight set, M1, 1 mg - 50 g, Knob, stainless steel fine turned, in wooden case

ENGLISH



 +49 7433 9933-0
 info@kern-sohn.com
 www.kern-sohn.com



FEATURES

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
353-420-210	Wood weight case 353-420-210

SERVICES

Model	Description
920-702	Marking (Last digits of serial number) KERN 920-702
962-601	DAkkS calibration certificate 962-601
962-601R	DAkkS Recalibration 962-601R

ENGLISH

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days