

Product brochure

07.08.2026

# KERN 343-054

Weight set, M1, 1 mg - 500 g, Knob, stainless steel fine turned, in plastic case


ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany


☎ +49 7433 9933-0  
✉ [info@kern-sohn.com](mailto:info@kern-sohn.com)  
🌐 [www.kern-sohn.com](http://www.kern-sohn.com)



FEATURES

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

PICTOGRAMS



ACCESSORIES

| Model       | Description                                                                                                             |
|-------------|-------------------------------------------------------------------------------------------------------------------------|
| 313-052-400 | Plastic carrying case 313-052-400 (for nominal values bis 500 g   for classes E1 - M3   for design Cylindrical/compact) |

SERVICES

| Model    | Description                                         |
|----------|-----------------------------------------------------|
| 920-705  | Marking (Last digits of serial number) KERN 920-705 |
| 962-604  | DAkKS calibration certificate 962-604               |
| 962-604R | DAkKS Recalibration 962-604R                        |

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days