

Product brochure

07.08.2026

# KERN 344-034

Weight set, M1, 1 g - 100 g, Knob, stainless steel fine turned, in plastic case

ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

+49 7433 9933-0  
info@kern-sohn.com  
www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

| Model       | Description                                                                                                             |
|-------------|-------------------------------------------------------------------------------------------------------------------------|
| 313-052-400 | Plastic carrying case 313-052-400 (for nominal values bis 500 g   for classes E1 - M3   for design Cylindrical/compact) |

SERVICES

| Model    | Description                                         |
|----------|-----------------------------------------------------|
| 920-711  | Marking (Last digits of serial number) KERN 920-711 |
| 962-616  | DAkKS calibration certificate 962-616               |
| 962-616R | DAkKS Recalibration 962-616R                        |

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days