

Product brochure

07.08.2026

KERN 344-07

Weight set, M1, 1 g - 2 kg, Knob, stainless steel fine turned, in wooden case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



PICTOGRAMS

Standard

Option



SERVICES

Model	Description
920-715	Marking (Last digits of serial number) KERN 920-715
962-620	DAkKS calibration certificate 962-620
962-620R	DAkKS Recalibration 962-620R

ENGLISH

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days