

# KERN 344-07

**KERN**

Weight set, M1, 1 g - 2 kg, Knob, stainless steelfine turned, in wooden case



## Optional cases

Case - Wood article number 353-470-200

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML M1
Product family	344

## Measuring System

Test weight - nominal value	1 g - 2 kg
OIML-compliant	✓
Weight OIML class	M1

## Construction

Dimension (W×D×H)	305×180×140 mm
Shape weight	Knob
Material	stainless steelfine turned
Number of weights	15
Case of weight set	in wooden case

## Services (optional)

Article number for conformity assessment	952-620
Article number for DAkkS calibration	962-620

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	305×180×140 mm
Shipping method	Parcel service
Net weight approx.	10 kg
Gross weight approx.	10 kg
Shipping weight	9,01 kg