

Product brochure

08.08.2026

KERN 344-076

Weight set, M1, 1 g - 2 kg, Knob, stainless steel fine turned, in aluminium case


ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
314-070-600	Aluminium weight case 314-070-600 (for nominal values 1 g - 2 kg for classes E1 - M2 for design Cylindrical/compact)

SERVICES

Model	Description
920-715	Marking (Last digits of serial number) KERN 920-715
962-620	DAkKS calibration certificate 962-620
962-620R	DAkKS Recalibration 962-620R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days