

Product brochure

08.08.2026

# KERN 344-086

Weight set, M1, 1 g - 5 kg, Knob, stainless steel fine turned, in aluminium case


ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

☎ +49 7433 9933-0  
✉ [info@kern-sohn.com](mailto:info@kern-sohn.com)  
🌐 [www.kern-sohn.com](http://www.kern-sohn.com)



PICTOGRAMS

Standard

Option

ACCESSORIES

| Model       | Description                                                                                                              |
|-------------|--------------------------------------------------------------------------------------------------------------------------|
| 314-080-600 | Aluminium weight case 314-080-600 (for nominal values 1 g - 5 kg   for classes E1 - M2   for design Cylindrical/compact) |

SERVICES

| Model    | Description                                         |
|----------|-----------------------------------------------------|
| 920-716  | Marking (Last digits of serial number) KERN 920-716 |
| 962-621  | DAkKS calibration certificate 962-621               |
| 962-621R | DAkKS Recalibration 962-621R                        |

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days