

KERN 344-086

KERN

Weight set, M1, 1 g - 5 kg, Knob, stainless steelfine turned, in aluminium case



Optional cases

Case - Aluminium article number 314-080-600

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML M1
Product family	344

Measuring System

Test weight - nominal value	1 g - 5 kg
OIML-compliant	✓
Weight OIML class	M1

Construction

Dimension (W×D×H)	460×260×200 mm
Shape weight	Knob
Material	stainless steelfine turned
Number of weights	16
Case of weight set	in aluminium case

Services (optional)

Article number for conformity assessment	952-621
Article number for DAkkS calibration	962-621

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	460×260×200 mm
Shipping method	Parcel service
Net weight approx.	13 kg
Gross weight approx.	18 kg
Shipping weight	18 kg