

Product brochure

07.08.2026

KERN 344-09

Weight set, M1, 1 g - 10 kg, Knob, stainless steel fine turned, in wooden case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
353-490-200	Wood weight case 353-490-200 (for nominal values 1 mg - 10 kg for classes F2 + M1 + M2 for design Cylindrical/compact)

SERVICES

Model	Description
920-717	Marking (Last digits of serial number) KERN 920-717
962-622	DAkKS calibration certificate 962-622
962-622R	DAkKS Recalibration 962-622R

KERN Pictogram Overview



Pallet shipment: The duration of internal product provision are two days



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days