

Product brochure

08.08.2026

KERN 344-096

Weight set, M1, 1 g - 10 kg, Knob, stainless steel fine turned, in aluminium case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
314-090-600	Aluminium weight case 314-090-600 (for nominal values 1 g - 10 kg for classes E1 - M2 for design Cylindrical/compact)

SERVICES

Model	Description
920-717	Marking (Last digits of serial number) KERN 920-717
962-622	DAkkS calibration certificate 962-622
962-622R	DAkkS Recalibration 962-622R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days