


Product brochure

08.08.2026

KERN 346-022-005

Weight container for rectangular weights, calibrated to OIML class M1

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

FEATURES

- Preconfigured weight container for testing high-load floor scales, pallet scales, pallet truck scales, crane scales, etc. This can also be used for storing the weights.
- This means the weight container and the weights can be placed on the balance in one go, saving you time and money
- The weight container is adjusted to OIML accuracy class M1. Other OIML accuracy classes are also available, please ask
- Weight of the weight container, OIML class M1: 20 kg
- Possible equipment, rectangular weights, OIML class M1: 5×20 kg
- Maximum total weight (weight container incl. weights): 120 kg
- Tolerance ± 1 g

PICTOGRAMS

Standard



ACCESSORIES

Model	Description
346-08	Test weight 346-08 (OIML Class M1, Nominal value 20 kg, Block, stainless steel)
346-78	Test weight 346-78 (OIML Class M1, Nominal value 20 kg, Block Eco, cast iron lacquered)
346-88	Test weight 346-88 (OIML Class M1, Nominal value 20 kg, Block, cast iron lacquered)
356-78	Test weight 356-78 (OIML Class M2, Nominal value 20 kg, Block Eco, cast iron lacquered)
356-88	Test weight 356-88 (OIML Class M2, Nominal value 20 kg, Block, cast iron lacquered)
366-78	Test weight 366-78 (OIML Class M3, Nominal value 20 kg, Block Eco, cast iron lacquered)
366-88	Test weight 366-88 (OIML Class M3, Nominal value 20 kg, Block, cast iron lacquered)

SERVICES

Model	Description
962-645	DAkkS calibration certificate 962-645
962-645R	DAkkS Recalibration 962-645R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day