


KERN 346-042-008

KERN

Weight container for rectangular weights, calibrated to OIML class M1



Pictograms

STANDARD



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Weight container

Measuring System

Test weight - nominal value	40 kg
OIML-compliant	✓
Weight OIML class	M1
Permissible tolerance according to OIML (+/-)	1,5 mg

Construction

Shape weight	weight container
Material	stainless steel

Services (optional)

Article number for DAkkS calibration	962-646
--------------------------------------	---------

Packing & Shipping

Delivery time	2 d
Dimensions packaging (W×D×H)	1000×500×500 mm
Net weight	40 kg
Shipping method	Freight
Net weight approx.	40 kg
Gross weight approx.	40 kg
Shipping weight	50 kg