


Product brochure

08.08.2026

KERN 346-065-009

Weight container for rectangular weights, calibrated to OIML class M1

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

FEATURES

- Preconfigured weight container for testing high-load floor scales, pallet scales, pallet truck scales, crane scales, etc. This can also be used for storing the weights.
- This means the weight container and the weights can be placed on the balance in one go, saving you time and money
- The weight container is adjusted to OIML accuracy class M1. Other OIML accuracy classes are also available, please ask
- Weight of the weight container, OIML class M1: 60 kg
- Possible equipment, rectangular weights, OIML class M1: 8×50 kg + 2×20 kg
- Maximum total weight (weight container incl. weights): 500 kg
- Tolerance ± 3 g

PICTOGRAMS

Standard



ACCESSORIES

Model	Description
346-08	Test weight 346-08 (OIML Class M1, Nominal value 20 kg, Block, stainless steel)
346-09	Test weight 346-09 (OIML Class M1, Nominal value 50 kg, Block, stainless steel)
346-78	Test weight 346-78 (OIML Class M1, Nominal value 20 kg, Block Eco, cast iron lacquered)
346-79	Test weight 346-79-2023a (OIML Class M1, Nominal value 50 kg, Block Eco, cast iron lacquered)
346-79	Test weight 346-79 (OIML Class M1, Nominal value 50 kg, Block Eco, cast iron lacquered)
346-88	Test weight 346-88 (OIML Class M1, Nominal value 20 kg, Block, cast iron lacquered)
346-89	Test weight 346-89 (OIML Class M1, Nominal value 50 kg, Block, cast iron lacquered)
356-78	Test weight 356-78 (OIML Class M2, Nominal value 20 kg, Block Eco, cast iron lacquered)
356-79	Test weight 356-79 (OIML Class M2, Nominal value 50 kg, Block Eco, cast iron lacquered)
356-88	Test weight 356-88 (OIML Class M2, Nominal value 20 kg, Block, cast iron lacquered)
356-89	Test weight 356-89 (OIML Class M2, Nominal value 50 kg, Block, cast iron lacquered)
366-78	Test weight 366-78 (OIML Class M3, Nominal value 20 kg, Block Eco, cast iron lacquered)
366-79	Test weight 366-79 (OIML Class M3, Nominal value 50 kg, Block Eco, cast iron lacquered)
366-88	Test weight 366-88 (OIML Class M3, Nominal value 20 kg, Block, cast iron lacquered)
366-89	Test weight 366-89 (OIML Class M3, Nominal value 50 kg, Block, cast iron lacquered)

SERVICES

Model	Description
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-646	DAkKS calibration certificate 962-646
962-646R	DAkKS Recalibration 962-646R

KERN Pictogram Overview



Pallet shipment: The duration of internal product provision are two days