

Product brochure

07.08.2026

KERN 346-08

Individual weight, OIML Class M1, 20 kg, Block, stainless steel

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PICTOGRAMS

Standard

Option





ACCESSORIES

Model	Description
317-280	Gloves 317-280 (cotton, 1 pair)
317-290	Gloves 317-290 (leather, cotton, 1 pair)
346-080-600	Aluminium weight case 346-080-600 (for nominal values 20 kg Block for classes F1 - M3 for design Block)

SERVICES

Model	Description
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-645	DAkkS calibration certificate 962-645
962-645R	DAkkS Recalibration 962-645R

ENGLISH

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days