


Product brochure

08.08.2026

KERN 347-022-005

Weight container for test weights, calibrated to OIML class M1

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

FEATURES

- Preconfigured weight container for testing high-load floor scales, pallet scales, pallet truck scales, crane scales, etc. This can also be used for storing the weights.
- This means the weight container and the weights can be placed on the balance in one go, saving you time and money
- The weight container is adjusted to OIML accuracy class M1. Other OIML accuracy classes are also available, please ask
- Weight of the weight container, OIML class M1: 20 kg
- Possible equipment, test weights, OIML class M1: 10×10 kg resp. 5×20 kg
- Maximum total weight (weight container incl. weights): 120 kg
- Tolerance ± 1 g

PICTOGRAMS

Standard



ACCESSORIES

Model	Description
327-141	Test weight 327-141 (OIML Class F1, Nominal value 10 kg, Check weight, stainless steel polished)
327-151	Test weight 327-151 (OIML Class F1, Nominal value 20 kg, Check weight, stainless steel polished)
337-141	Test weight 337-141 (OIML Class F2, Nominal value 10 kg, Check weight, stainless steel fine turned)
337-151	Test weight 337-151 (OIML Class F2, Nominal value 20 kg, Check weight, stainless steel fine turned)
347-141	Test weight 347-141 (OIML Class M1, Nominal value 10 kg, Check weight, stainless steel fine turned)
347-151	Test weight 347-151 (OIML Class M1, Nominal value 20 kg, Check weight, stainless steel fine turned)

SERVICES

Model	Description
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-645	DAkkS calibration certificate 962-645
962-645R	DAkkS Recalibration 962-645R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day