

Product brochure

07.08.2026

KERN 347-05

Individual weight, OIML Class M1, 20 g, Knob, stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
315-250	Tweezers 315-250 (stainless steel 150 mm)
317-050-600	Aluminium weight case 317-050-600 (for nominal values 20 g for classes E1 - M3 for design Cylindrical/compact)
337-050-200	Wood weight case 337-050-200 (for nominal values 20 g for classes F2 + M1 for design Cylindrical/compact)
347-050-400	Plastic weight case 347-050-400 (for nominal values 10 g - 20 g for classes F1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-635	DAkkS calibration certificate 962-635
962-635R	DAkkS Recalibration 962-635R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days