

# KERN 347-050-400

**KERN**

Plastic box for individual weights (Cylindrical/compact, 10 g - 20 g)



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Plastic box

## Measuring System

Case suitable for OIML classes	F1 - M3
Case suitable for nominal weight values	10 g - 20 g
Suitable weight shape for case	Cylindrical/compact

## Construction

Dimension (Ø×H)	24×51 mm
Material	plastic
Case - Smallest possible weight	20 g
Case - Highest possible weight	20 g

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	24×24×51 mm
Net weight	0,02 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

## Pictograms

STANDARD

