


KERN 347-062-022

KERN

Weight container for test weights, calibrated to OIML class M1



Pictograms

STANDARD



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Weight container
Product family	347

Measuring System

Test weight - nominal value	60 kg
OIML-compliant	✓
Weight OIML class	M1
Permissible tolerance according to OIML (+/-)	3 g

Construction

Shape weight	Weight container
Material	stainless steel

Services (optional)

Article number for DAkkS calibration	962-646
--------------------------------------	---------

Packing & Shipping

Delivery time	2 d
Dimensions packaging (W×D×H)	1000×500×500 mm
Shipping method	Freight
Net weight approx.	60 kg
Gross weight approx.	60 kg
Shipping weight	60 kg