

KERN 347-07

Individual weight, OIML Class M1, 100 g, Knob, stainless steelfine turned



Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	22×22×39 mm
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg

Optional cases

Case - Plastic article number	347-070-400
Case - Aluminium article number	317-070-600
Case - Wood article number	337-070-200

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

Measuring System

Conformity assessment possible	✓
Test weight - volume (cm ³)	12,5
Test weight - volume uncertainty (cm ³)	0,1563
Test weight - Density (kg/m ³)	7950
Test weight - Density uncertainty (kg/m ³)	140
Test weight - Density [Min] (OIML)	2600
Test weight - nominal value	100 g
OIML-compliant	✓
Weight OIML class	M1
Permissible tolerance according to OIML (+/-)	5 mg

Construction

Shape weight	Knob
Dimensions (Ø×H)	22×39 mm
Material	stainless steelfine turned

Services (optional)

Article number for conformity assessment	952-637
Article number for DAkkS calibration	962-637