

Product brochure

07.08.2026

KERN 347-14

Individual weight, OIML Class M1, 10 kg, Knob, stainless steel fine turned

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
317-140-600	Aluminium weight case 317-140-600 (for nominal values 10 kg for classes E1 - M3 for design Cylindrical/compact)
317-280	Gloves 317-280 (cotton, 1 pair)
317-290	Gloves 317-290 (leather, cotton, 1 pair)
337-140-200	Wood weight case 337-140-200 (for nominal values 10 kg for classes F2 + M1 for design Cylindrical/compact)
347-140-400	Plastic weight case 347-140-400 (for nominal values 10 kg for classes F1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
962-644	DAkKS calibration certificate 962-644
962-644R	DAkKS Recalibration 962-644R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days