

Product brochure

07.08.2026

# KERN 347-161

Individual weight, OIML Class M1, 50 kg, Check weight, stainless steel fine turned

ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

+49 7433 9933-0  
info@kern-sohn.com  
www.kern-sohn.com

FEATURES

- stackable

PICTOGRAMS

Standard

Option





ACCESSORIES

| Model       | Description                                                                                             |
|-------------|---------------------------------------------------------------------------------------------------------|
| 317-290     | Gloves 317-290 (leather, cotton, 1 pair)                                                                |
| 337-161-200 | Wood weight case 337-161-200 (for nominal values 50 kg   for classes F2 + M1   for design Check weight) |

SERVICES

| Model    | Description                                                         |
|----------|---------------------------------------------------------------------|
| 920-350  | Marking by laser engraving (individual marking) KERN 920-350        |
| 920-351  | Marking of test weights (last digits of serial number) KERN 920-351 |
| 962-646  | DAkKS calibration certificate 962-646                               |
| 962-646R | DAkKS Recalibration 962-646R                                        |

## KERN Pictogram Overview



Pallet shipment: The duration of internal product provision are two days



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days