

# KERN 353-420-200

**KERN**

Wooden box for weight sets (Cylindrical/compact, 1 mg - 50 g)



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

## Measuring System

Case suitable for OIML classes	F2 + M1 + M2
Case suitable for nominal weight values	1 mg - 50 g
Suitable weight shape for case	Cylindrical/compact

## Construction

Dimension (W×D×H)	129×86×40 mm
Material	beech
Case - Smallest possible weight	1 mg
Case - Highest possible weight	50 g

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	129×86×40 mm
Net weight	0,2 kg
Shipping method	Parcel service
Net weight approx.	0,20 kg
Gross weight approx.	0,20 kg
Shipping weight	0,2 kg

## Pictograms

STANDARD

