

KERN 353-430-200

KERN

Wooden box for weight sets (Cylindrical/compact, 1 mg - 100 g)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for OIML classes	F2 + M1 + M2
Case suitable for nominal weight values	1 mg - 100 g
Suitable weight shape for case	Cylindrical/compact

Construction

Dimension (W×D×H)	151×91×49 mm
Material	beech
Case - Smallest possible weight	1 mg
Case - Highest possible weight	100 g

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	151×91×49 mm
Net weight	0,25 kg
Shipping method	Parcel service
Net weight approx.	0,25 kg
Gross weight approx.	0,25 kg
Shipping weight	0,25 kg

Pictograms

STANDARD

