

Product brochure

07.08.2026

KERN 353-440-200

Wooden box for weight sets (Cylindrical/compact, 1 mg - 200 g)

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

PICTOGRAMS

Standard



ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
318-270	Brush 318-270
353-420-210	Wood weight case 353-420-210

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350

ENGLISH

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day