

KERN 353-450-200

Wooden box for weight sets (Cylindrical/compact, 1 mg - 500 g)



Category

| | |
|------------------|---------------|
| Brand | KERN |
| Product category | Cases & Boxes |
| Product group | Wooden box |

Measuring System

| | |
|---|---------------------|
| Case suitable for OIML classes | F2 + M1 + M2 |
| Case suitable for nominal weight values | 1 mg - 500 g |
| Suitable weight shape for case | Cylindrical/compact |

Construction

| | |
|---------------------------------|--------------|
| Dimension (W×D×H) | 168×98×78 mm |
| Material | beech |
| Case - Smallest possible weight | 1 mg |
| Case - Highest possible weight | 500 g |

Packing & Shipping

| | |
|------------------------------|----------------|
| Delivery time | 1 d |
| Dimensions packaging (W×D×H) | 168×98×78 mm |
| Net weight | 0,36 kg |
| Shipping method | Parcel service |
| Net weight approx. | 0,40 kg |
| Gross weight approx. | 0,40 kg |
| Shipping weight | 0,36 kg |

Pictograms

STANDARD

