

KERN 353-470-200

KERN

Wooden box for weight sets (Cylindrical/compact, 1 mg - 2 kg)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for OIML classes	F2 + M1 + M2
Case suitable for nominal weight values	1 mg - 2 kg
Suitable weight shape for case	Cylindrical/compact

Construction

Dimension (W×D×H)	305×180×140 mm
Material	beech
Case - Smallest possible weight	1 mg
Case - Highest possible weight	2 kg

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	305×180×140 mm
Net weight	3,2 kg
Shipping method	Parcel service
Net weight approx.	3,2 kg
Gross weight approx.	3,6 kg
Shipping weight	3,5 kg

Pictograms

STANDARD

