

# KERN 354-046

# KERN

Weight set, M2, 1 g - 200 g, Knob, stainless steel fine turned, in aluminium case



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML M2
Product family	354

## Measuring System

Test weight - nominal value	1 g - 200 g
OIML-compliant	✓
Weight OIML class	M2

## Construction

Dimension (W×D×H)	190×130×80 mm
Shape weight	Knob
Material	stainless steel fine turned
Case of weight set	in aluminium case

## Services (optional)

Article number for DAkkS calibration	962-617
--------------------------------------	---------

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	190×130×80 mm
Shipping method	Parcel service
Net weight approx.	0,90 kg
Gross weight approx.	1,2 kg
Shipping weight	1,1 kg

## Optional cases

Case - Aluminium article number	314-040-600
---------------------------------	-------------