

KERN 354-056

KERN

Weight set, M2, 1 g - 500 g, Knob, stainless steel fine turned, in aluminium case



Optional cases

Case - Aluminium article number 314-050-600

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML M2
Product family	354

Measuring System

Test weight - nominal value	1 g - 500 g
OIML-compliant	✓
Weight OIML class	M2

Construction

Dimension (W×D×H)	210×140×98 mm
Shape weight	Knob
Material	stainless steel fine turned
Case of weight set	in aluminium case

Services (optional)

Article number for DAkkS calibration	962-618
--------------------------------------	---------

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	210×140×98 mm
Net weight	1,6 kg
Shipping method	Parcel service
Net weight approx.	1,6 kg
Gross weight approx.	2,0 kg
Shipping weight	2 kg