

KERN 354-420-200

KERN

Wooden box for weight sets (Cylindrical/compact, 1 g - 50 g)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for OIML classes	F2 + M1 + M2 + M3
Case suitable for nominal weight values	1 g - 50 g
Suitable weight shape for case	Cylindrical/compact

Construction

Dimension (W×D×H)	121×53×40 mm
Material	beech
Case - Smallest possible weight	1 g
Case - Highest possible weight	50 g

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	120×50×40 mm
Net weight	0,15 kg
Shipping method	Parcel service
Net weight approx.	0,15 kg
Gross weight approx.	0,15 kg
Shipping weight	0,15 kg

Pictograms

STANDARD

