

Product brochure

07.08.2026

KERN 354-430-200

Wooden box for weight sets (Cylindrical/compact, 1 g - 100 g)

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PICTOGRAMS

Standard



ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350

ENGLISH

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day