

# KERN 354-450-200

Wooden box for weight sets (Cylindrical/compact, 1 g - 500 g)



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

## Measuring System

Case suitable for OIML classes	F2 + M1 + M2 + M3
Case suitable for nominal weight values	1 g - 500 g
Suitable weight shape for case	Cylindrical/compact

## Construction

Dimension (W×D×H)	183×82×80 mm
Material	beech
Case - Smallest possible weight	1 g
Case - Highest possible weight	500 g

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	170×80×80 mm
Net weight	0,61 kg
Shipping method	Parcel service
Net weight approx.	0,65 kg
Gross weight approx.	0,65 kg
Shipping weight	0,61 kg

## Pictograms

### STANDARD

