

KERN 362-97

KERN

Weight set, M3, 1 g - 2 kg, stainless steel (1 g - 50 g) and cast iron lacquered (≥ 100 g), button (≤ 50 g, ≥ 500 g) and cylinders (100 g, 200 g), in wood block



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML M3
Product family	362

Measuring System

Conformity assessment possible	no
Test weight - nominal value	1 g - 2 kg
OIML-compliant	✓
Weight OIML class	M3

Construction

Dimension (W×D×H)	340×120×60 mm
Shape weight	Knob
Material	cast iron lacquered
Number of weights	15
Case of weight set	in wood block

Services (optional)

Article number for DAkkS calibration	962-620
--------------------------------------	---------

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	340×120×60 mm
Shipping method	Parcel service
Net weight approx.	8 kg
Gross weight approx.	8 kg
Shipping weight	7,3 kg

Optional cases

Case - Wood article number	362-970-300
----------------------------	-------------