

# KERN 362-980-300

**KERN**

Wooden box for weight sets (Cylindrical/compact, 1 g - 5 kg)



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

## Measuring System

Case suitable for OIML classes	M3
Case suitable for nominal weight values	1 g - 5 kg
Suitable weight shape for case	Cylindrical/compact

## Construction

Dimension (W×D×H)	390×150×70 mm
Material	beech
Case - Smallest possible weight	1 g
Case - Highest possible weight	5 kg

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	390×150×70 mm
Net weight	2,2 kg
Shipping method	Parcel service
Net weight approx.	2,2 kg
Gross weight approx.	2,2 kg
Shipping weight	2,2 kg

## Pictograms

### STANDARD

